AD SYRINGE WITH NEEDLE

Product Type

AD syringe with needle	Needle Size	Basic UDI-DI: 611800222ADSFXV9
Luer 3 parts Syringe 0.5ml	23G 24G 25G	

Product Components

Name	Raw material
Barrel	PP - Polypropylene
Plunger	PP - Polypropylene
Gasket	Thermoplastic Elastomer, Latex free
Hypodermic Needle	Stainless Steel, SS304, PP - Polypropylene
Packaging	Medical grade paper and transparent film for EtO (individually blister packed)

Product Description

- The auto-disable (AD) syringe is a sterile, single-use medical device designed for immunization and medication administration. It consists of a calibrated plastic barrel with a smoothly moving plunger, thanks to an integrated gasket, and a fixed, capped hypodermic needle for user protection.
- Made from plastic and thermoplastic elastomer materials, it features a built-in re-use prevention mechanism that activates automatically at the start of injection, locking or breaking the plunger after a single use, eliminating the risk of reuse.
- The AD syringe ensures patient safety by preventing contamination and cross-infection, making it a crucial tool in vaccination programs and medical treatments.

Intended Use

- AD Syringe with needle is used for immunization with a fixed dose. It is used for administering medications/Vaccines. It consists of a transparent barrel to ensure controlled administration of medicines with graduation, a plunger with seal/gasket and a fixed hypodermic needle (capped for user protection).
- Auto disable syringe is a sterile (fluid path), Non-Pyrogenic & Non-toxic product. This device is used for medical purpose injection and aspiration of fluids from vials, ampoules, and parts of body below the surface of the skin.

Intended User

The medical device shall be used by Healthcare Professional.

Patient Population

AD Syringe for fixed-dose immunization is designed for use with individuals of all ages, from newborns to adults and the elderly.


AD SYRINGE WITH NEEDLE

Instructions for Use


DON'T PUSH THE PLUNGER FORWARD BEFORE FILLING THE MEDICATION.

Step 1: Open the pack and take out the syringe.
Do not push the plunger forward. Otherwise, the syringe will get locked.


Step 2: Insert the needle into the vial and withdraw the fluid/medication.

Step 3: Flick the barrel to move any air bubbles to the top.
Push the plunger to the desired dosage and expel the air bubbles.

Step 4: Inject it to the patient until the liquid is completely aspirated.

Step 5: Continue with the normal injection procedure, including aspiration to ensure the vein is not punctured. If the plunger is pulled back after the injection, it will break.

Step 6: Dispose of the syringe in a sharps container.

The electronic version of "IFU/ADSFX/001 - Instructions for use" is available at www.immedicare.ma.

AD SYRINGE WITH NEEDLE

Compatibility with other medical devices

Suitable for use with connectors "EN ISO 80369-7 Small-bore connectors for liquids and gases in healthcare applications - Part 7: Connectors for intravenous or hypodermic applications".

Indications

- AD Syringe with needle is intended for the injection of intravenous and/or intramuscular fluids immediately after filling.
- These devices play a role in the treatment and alleviation of diseases by administering fixed dosages of medicinal fluids into the body.

Contraindications

- The device must be used only as indicated in its intended use.
- It should not be used for the administration of blood and blood components.
- Additionally, the device is not suitable for large volume fluids and fluid derivatives, nor for high-viscosity fluids.

Cautions

- Do not use the product if the package is opened and damaged.
- Visually inspect that the content of the package is complete, and the package is undamaged.
- Do not use the product after the expiration date printed on the primary packaging.
- After injection or aspiration, dispose of the product as medical waste.
- The product is for single use only. It cannot be reused or re-sterilized.
- Do not touch the needle tip to prevent contamination/needle sticks. Skin puncture with a contaminated needle may result in serious illnesses such as hepatitis, AIDS and known/unknown diseases.
- If the sterile packaging is not used immediately after opening, there is a risk of contamination. Product sterility cannot be maintained, do not use the product.
- Reuse of the product may cause infections, cross-contamination and sepsis.
- Sterile syringes are designed to be used immediately after being filled, so do not store the medicine for a long time.
- If the wrapping or cap on the end of the syringe is broken or missing, discard it in your sharps container and use a new needle and syringe.

Disclaimer

- The Medical device is supplied with limited functionality, I'M MEDICARE does not warrant any foreseeable risk resulting from improper/ misuse of the medical device
- Device should not be used out of its instructions for use. In case the device is used out of its instructions described, the manufacturer/supplier must be informed for any usage outside of its recommendations.
- If any serious incident occurs, it should be reported to the manufacturer and the competent authority.

For any feedback or inquiries, customers/users can contact us at:

Email: contact@immedicare.ma

Customer Care No: +212 524 32 08 77

AD SYRINGE WITH NEEDLE

Adverse Event

- The device is supplied with limited functionality, and there are no foreseeable risks associated with its use.
- Improper use may cause redness, swelling, bruising, or other adverse effects.

Device Malfunction/Failure

In case of a device malfunction or failure, follow these steps:

- Discontinue use immediately if you notice any leakage, difficulty in activation, or failure to perform as intended.
- Do not attempt to repair the device yourself.
- Dispose of the device safely following the provided disposal guidelines.
- Report the malfunction to the manufacturer or an authorized representative for further investigation.

Storage Conditions


- Store in a cool and dry place.
- Keep it away from direct sunlight and moist environment.

The shelf life of product is 5 years.

AD SYRINGE WITH NEEDLE

Description of Symbols Used

Symbol	Description	Symbol	Description
	Manufactured by		Keep away from sunlight
	Date of manufacturing		Keep dry
	Use-by / Expiry date		Temperature limitation
	Lot number / Batch number		Moisture limitation
	Product reference		Do not use if the package is damaged
	Medical device		Single sterile barrier system
	Unique device identifier		Do not reuse
	European Authorized Representative		Latex free
	CE marking		Non-pyrogenic
	Sterilized by Ethylene oxide		Consult for IFU
	Do not re-sterilize		Do not use hook
	Caution / Warning (Read IFU before use)		Stock and store in upright condition
	Fragile. Handle with care		Re-use prevention

IM MEDICARE

Address: Lot 21, Parc Industriel Sidi Bou Othmane, 43303, Morocco

Email: contact@immedicare.ma

Phone number: +212 524 32 08 77

Fax: +212 524 32 06 29

Website: www.immedicare.ma

Name: EUCEREP

Address: Roald Dahllaan 33, 5629MC, Eindhoven, The Netherlands

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