
HYPODERMIC NEEDLE

Product Type

Hypodermic Needle								Basic UDI-DI: 611800222HYPONDLEL
16G	17G	18G	20G	21G	22G	23G	24G	
25G	26G	28G	29G	30G	31G	32G		

Product Components

Name	Raw material
Hypodermic Needle	Stainless Steel, SS304, PP – Polypropylene
Protective/Safety Cap	PP - Polypropylene
Packing Material	Medical grade paper and transparent film for EtO (individually blister packed)

Product Description

- A non-dedicated, sharp bevel-edged, hollow tubular device designed to inject substances into the body or extract fluids from the body while connected to a noninvasive device. Some variants include a protective shielding mechanism (e.g., safety cap) to prevent sharp injuries.
- It is not specifically intended for blood collection or arterial puncture and lacks fixation wings, distinguishing it from blood collection, intra-arterial, or venous butterfly needles.
- This is a single-use device, ensuring safety and reducing the risk of contamination.

Intended Use

- Hypodermic needles are used to inject substances into the body or extract fluids from the body while connected to a noninvasive device (e.g., syringe, secondary medication set).
- A hypodermic needle offers a fast, safe & easy option for rapid delivery of liquids to the patient when the injected substance cannot be ingested orally, as it would not be absorbed or would harm organs.

Intended User

The medical device shall be used by Healthcare Professional.

Patient Population

Hypodermic needles are designed for use with individuals of all ages, from newborns to adults and the elderly.

Instructions for Use

- Open the packaging by stripping it from the opening place indicated by the arrow sign.
- Visually inspect that the content of the package is complete, and the package is undamaged.
- The product packaging consists of a sterile needle.

HYPODERMIC NEEDLE

- Without removing the protective cap of the needle, firmly secure the syringe and needle together, pushing the needle toward the syringe tip and turning it clockwise.
- Carefully remove the protective cap of the needle without damaging the needle. For safety hypodermic needles, carefully remove the safety cap of the needle without damaging the needle.
- Inject the liquid to be injected into the intravenously, intramuscularly, subcutaneously or intradermally, according to the application procedure.

The electronic version of "IFU/HYPONDL/001 - Instructions for use" is available at www.immedicare.ma.

Compatibility with other medical devices

Suitable for use with connectors "EN ISO 80369-7 Small-bore connectors for liquids and gases in healthcare applications - Part 7: Connectors for intravenous or hypodermic applications".

Indications

- Hypodermic needles made of empty, sterile and single use stainless steel are used by end users for the injection or withdrawal of both body fluids (blood, etc.) and medical fluids such as medicines.
- It may also be used for liquid/drug preparation.

Contraindications

- Hypodermic needles are not used in more than one patient. The product can be used in a one-time sterile condition with a single patient.
- Cleaning or re-sterilization cannot be performed. It is not intended for use outside the specified intended use. After usage, dispose of the product as medical waste.

Cautions

- Do not use the product if the package is opened and damaged.
- Visually inspect that the content of the package is complete, and the package is undamaged.
- Do not use the product after the expiration date printed on the primary packaging.
- Before use, check that the hypodermic needle is not blocked, and the flow is not obstructed.
- For safety hypodermic needles, activate it away from yourself and others with the one-hand technique. For a high level of safety, use only the finger pad area on the large rough tissue when activating the safety cap. If liquid remains on the needle during the activation of the safety cap or after the injection, it may cause a small amount of splash.
- For safety hypodermic needles, to help avoid HIV (AIDS), HBV (Hepatitis) and other infectious diseases due to accidental needlesticks, activate the protective mechanism immediately after use.
- After injection or aspiration, dispose of the product as medical waste.
- The product is for single use only. It cannot be reused, reworked or re-sterilized.
- If the sterile packaging is not used immediately after opening, there is a risk of contamination. Product sterility cannot be maintained, do not use the product.
- Reuse of the product may cause infections, cross-contamination and sepsis.

HYPODERMIC NEEDLE

Disclaimer

- The Medical device is supplied with limited functionality, I'M MEDICARE does not warrant any foreseeable risk resulting from improper/ misuse of the medical device
- Device should not be used out of its instructions for use. In case the device is used out of its instructions described, the manufacturer/supplier must be informed for any usage outside of its recommendations.
- If any serious incident occurs, it should be reported to the manufacturer and the competent authority.

For any feedback or inquiries, customers/users can contact us at:

Email: contact@immedicare.ma

Customer Care No: +212 524 32 08 77

Adverse Event

- The device is supplied with limited functionality, and there are no foreseeable risks associated with its use.
- Improper use may cause redness, swelling, bruising, or other adverse effects.

Device Malfunction/Failure

In case of a device malfunction or failure, follow these steps:

- Discontinue use immediately if you notice any leakage, difficulty in activation, or failure to perform as intended.
- Do not attempt to repair the device yourself.
- Dispose of the device safely following the provided disposal guidelines.
- Report the malfunction to the manufacturer or an authorized representative for further investigation.



Storage Conditions

- Store in a cool and dry place.
- Keep it away from direct sunlight and moist environment.
- The shelf life of product is 5 years.

HYPODERMIC NEEDLE

Description of Symbols Used

Symbol	Description	Symbol	Description
	Manufactured by		Keep away from sunlight
	Date of manufacturing		Keep dry
	Use-by / Expiry date		Temperature limitation
	Lot number / Batch number		Moisture limitation
	Product reference		Do not use if the package is damaged
	Medical device		Single sterile barrier system
	Unique device identifier		Do not reuse
	European Authorized Representative		Latex free
	CE marking		Non-pyrogenic
	Sterilized by Ethylene oxide		Consult for IFU
	Do not re-sterilize		Do not use hook
	Caution / Warning (Read IFU before use)		Stock and store in upright condition
	Fragile. Handle with care		

I'M MEDICARE

Address: Lot 21, Parc Industriel Sidi Bou Othmane, 43303, Morocco

Email: contact@immedicare.ma

Phone number: +212 524 32 08 77

Fax: +212 524 32 06 29

Website: www.immedicare.ma

Name: EUCEREP


Address: Roald Dahllaan 33, 5629MC, Eindhoven, The Netherlands

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