STERILE HYPODERMIC SYRINGE (WITH/WITHOUT NEEDLE)

Product Type

SOLOJEKT®	Needle Size	Basic UDI-DI: 611800222HYPOSRGD
Luer 3 parts Syringe 1 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50/60 ml	18G 21G 22G 23G 24G 25G	

Product Components

Name	Raw material
Barrel	PP - Polypropylene
Plunger	PP - Polypropylene
Gasket	Thermoplastic Elastomer, Latex free
Hypodermic Needle	Stainless Steel, SS304, PP - Polypropylene
Packaging	Medical grade paper and transparent film for EtO (individually blister packed)

Product Description

- A sterile device consisting of a calibrated barrel (cylinder) with plunger (attached with a gasket allowing smooth plunger movement) intended for the aspiration of fluids or the injection immediately after filling, a non-retractable needle is included (optional), either attached or detached (capped for user protection).
- The syringe is typically made of plastic and Thermoplastic Elastomer materials
- This is a single-use device.

Intended Use

Sterile Hypodermic Syringe for single use intended for the aspiration of fluids or the injection immediately after filling, it is mostly used for the administration of medications or for blood sampling.

Intended User

The medical device shall be used by Healthcare Professional.

Patient Population

Sterile hypodermic syringe is designed for use with individuals of all ages, from newborns to adults and the elderly.

Instructions for Use

- Inspect package for integrity and expiry and then remove the syringe from package using aseptic technique.
- Select the proper size of syringe & needle for injection.
- Carefully select and clean the venipuncture or muscle site.
- Be careful to pull the needle cap straight off to avoid the needle point damage.

STERILE HYPODERMIC SYRINGE (WITH/WITHOUT NEEDLE)

- Push and twist the needle clockwise on the syringe tip to seal it firmly and ensure a tight fitting.
- Draw the prescribed dosage into the barrel.
- Push the plunger carefully to take out any air bubbles.
- Inject the patient carefully with the needle.
- Push the Plunger thoroughly to the end of complete injection.
- Withdraw the needle from patient.
- Do not recap the needle.
- Pull the plunger to full extended position and break it to prevent anyone from using it again.
- Discard the needle after single use to prevent anyone from using it again

The electronic version of "IFU/HYPOSRG/001 - Instructions for use" is available at www.immedicare.ma.

Compatibility with other medical devices

Suitable for use with connectors "EN ISO 80369-7 Small-bore connectors for liquids and gases in healthcare applications - Part 7: Connectors for intravenous or hypodermic applications".

Indications

- Hypodermic syringes with needle are devices used to administer medical products and liquids.
- Syringes are designed to be used for intravenous, intramuscular, subcutaneous and taking blood samples.

Contraindications

- Hypodermic syringes with needle should not be used in more than one patient. The product can be used in a one-time sterile condition with a single patient.
- Cleaning or re-sterilization is not possible. It is not intended for use other than its stated intended use.
- Dispose of as medical waste after use.

Cautions

- Do not use the product if the package is damaged.
- Visually inspect that the content of the package is complete, and the package is undamaged.
- Do not use the product after the expiration date printed on the primary packaging.
- After injection or aspiration, dispose of the product as medical waste.
- The product is for single use only. It cannot be reused or re-sterilized.
- Do not touch the needle tip to prevent contamination/ needle sticks. Skin puncture with a contaminated needle may result in serious illnesses such as hepatitis, AIDS and known/unknown diseases.
- If the sterile packaging is not used immediately after opening, there is a risk of contamination. Product sterility cannot be maintained, do not use the product.
- Reuse of the product may cause infections, cross-contamination and sepsis.
- Sterile syringes are designed to be used immediately after being filled, so do not store the medicine for a long time.

STERILE HYPODERMIC SYRINGE (WITH/WITHOUT NEEDLE)

Disclaimer

- The Medical device is supplied with limited functionality, I'M MEDICARE does not warrant any foreseeable risk resulting from improper/ misuse of the medical device
- Device should not be used out of its instructions for use. In case the device is used out of its instructions described, the manufacturer/supplier must be informed for any usage outside of its recommendations.
- If any serious incident occurs, it should be reported to the manufacturer and the competent authority.

For any feedback or inquiries, customers/users can contact us at:

Email: contact@immedicare.ma

Customer Care No: +212 524 32 08 77

Adverse Event

- The device is supplied with limited functionality, and there are no foreseeable risks associated with its use.
- Improper use may cause redness, swelling, bruising, or other adverse effects.

Device Malfunction/Failure

In case of a device malfunction or failure, follow these steps:

- Discontinue use immediately if you notice any leakage, difficulty in activation, or failure to perform as intended.
- Do not attempt to repair the device yourself.
- Dispose of the device safely following the provided disposal guidelines.
- Report the malfunction to the manufacturer or an authorized representative for further investigation.

Storage Conditions

- Store in a cool and dry place.
- Keep it away from direct sunlight and moist environment.
- The shelf life of product is 5 years.

STERILE HYPODERMIC SYRINGE (WITH/WITHOUT NEEDLE)

Description of Symbols Used

Symbol	Description	Symbol	Description
	Manufactured by		Keep away from sunlight
	Date of manufacturing		Keep dry
	Use-by / Expiry date		Temperature limitation
	Lot number / Batch number		Moisture limitation
	Product reference		Do not use if the package is damaged
	Medical device		Single sterile barrier system
	Unique device identifier		Do not reuse
	European Authorized Representative		Latex free
	CE marking		Non-pyrogenic
	Sterilized by Ethylene oxide		Consult for IFU
	Do not re-sterilize		Do not use hook
	Caution / Warning (Read IFU before use)		Stock and store in upright condition
	Fragile. Handle with care		

I'M MEDICARE

Address: Lot 21, Parc Industriel Sidi Bou Othmane, 43303, Morocco

Email: contact@immedicare.ma

Phone number: +212 524 32 08 77

Fax: +212 524 32 06 29

Website: www.immedicare.ma

Name: EUCEREP


Address: Roald Dahllaan 33, 5629MC, Eindhoven, The Netherlands

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