INSULIN SYRINGE

Product Type

SOLOJEKT®	Needle Size	Basic UDI-DI: 611800222SYRINSBL
Luer 3 parts Syringe 0.5 ml, 1 ml	29G 30G 31G	

Product Components

Name	Raw material
Barrel	PP - Polypropylene
Plunger	PP - Polypropylene
Gasket	Thermoplastic Elastomer, Latex free
Hypodermic Needle	Stainless Steel, SS304, PP - Polypropylene
Packaging	Medical grade paper and transparent film for EtO (individually blister packed)

Product Description

Insulin syringes are designed for administering insulin injections to patients subcutaneously. Each insulin syringe consists of a graduated barrel, a plunger with a gasket, and a needle/hub assembly. The syringes are sterile, single-use, and non-pyrogenic. This device operates on principles similar to those of piston syringes.

Intended Use

A sterile device intended to administer insulin injection to a patient subcutaneously. It can be used by healthcare personnel or patients. It will allow for the aspiration of the pharmaceutical from its container for direct administration to the patient subcutaneously. This is a single-use device.

Intended User

Device should be used by a Doctor/Registered Practitioner or Paramedic staff or a patient, as guided by healthcare professional

Patient Population

The single uses Syringes for insulin are designed for administration of insulin to all individuals (newborn to adults to old age). The insulin syringe may be used by a patient or by a healthcare professional, either in the home setting, in a hospital, or a clinical setting.

INSULIN SYRINGE

Instructions for Use

Step 1: Open package by peeling back paper and remove syringe carefully.

Step 2: Remove protective needle sheath/cap carefully to avoid needle point damage.

Step 3: Draw up medication and administer injection using an aseptic technique in accordance with your facility's protocol.

Step 4: Once the insulin syringe is used, immediately dispose of syringe in an approved sharps disposal container.

The electronic version of "IFU/SYRINS/001 - Instructions for use" is available at www.immedicare.ma.

Compatibility with other medical devices

Suitable for use with connectors "EN ISO 80369-7 Small-bore connectors for liquids and gases in healthcare applications - Part 7: Connectors for intravenous or hypodermic applications".

Indications

Insulin syringes with needles are devices used to administer medical products and liquids.

Contraindications

- Device must not be used other than as indicated in the intended use.
- Device should not be used for the administration of blood and blood components.
- Device should not be used for large volume fluids and fluid derivatives.
- Device should not be used for the administration of high viscous fluids.
- Device should not be used on patients with known hypersensitivity to device material.

INSULIN SYRINGE

Cautions

- The device is for single use only.
- The device must not be reused. Reuse of the device may cause serious infection, cross-contamination, and transmission of bacteria, viruses, and blood-related diseases.
- Ensure before use that the device/packaging is not tampered with or damaged.
- The device is for manual use only.

Disclaimer

- The Medical device is supplied with limited functionality, I'M MEDICARE does not warrant any foreseeable risk resulting from improper/ misuse of the medical device
- Device should not be used out of its instructions for use. In case the device is used out of its instructions described, the manufacturer/supplier must be informed for any usage outside of its recommendations.
- If any serious incident occurs, it should be reported to the manufacturer and the competent authority.

For any feedback or inquiries, customers/users can contact us at:

Email: contact@immedicare.ma

Customer Care No: +212 524 32 08 77

Adverse Event

- The device is supplied with limited functionality, and there are no foreseeable risks associated with its use.
- Improper use may cause redness, swelling, bruising, or other adverse effects.

Device Malfunction/Failure

In case of a device malfunction or failure, follow these steps:

- Discontinue use immediately if you notice any leakage, difficulty in activation, or failure to perform as intended.
- Do not attempt to repair the device yourself.
- Dispose of the device safely following the provided disposal guidelines.
- Report the malfunction to the manufacturer or an authorized representative for further investigation.


Storage Conditions

- Store in a cool and dry place.
- Keep it away from direct sunlight and moist environment.
- The shelf life of product is 5 years.

INSULIN SYRINGE

Description of Symbols Used

Symbol	Description	Symbol	Description
	Manufactured by		Keep away from sunlight
	Date of manufacturing		Keep dry
	Use-by / Expiry date		Temperature limitation
	Lot number / Batch number		Moisture limitation
	Product reference		Do not use if the package is damaged
	Medical device		Single sterile barrier system
	Unique device identifier		Do not reuse
	European Authorized Representative		Latex free
	CE marking		Non-pyrogenic
	Sterilized by Ethylene oxide		Consult for IFU
	Do not re-sterilize		Do not use hook
	Caution / Warning (Read IFU before use)		Stock and store in upright condition
	Fragile. Handle with care		

I'M MEDICARE

Address: Lot 21, Parc Industriel Sidi Bou Othmane, 43303, Morocco

Email: contact@immedicare.ma

Phone number: +212 524 32 08 77

Fax: +212 524 32 06 29

Website: www.immedicare.ma

Name: EUCEREP


Address: Roald Dahllaan 33, 5629MC, Eindhoven, The Netherlands

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