
HYPODERMIC AD/RUP SYRINGE (WITH/WITHOUT NEEDLE)

Product Type

SOLOJEKT AD®	Needle Size	Basic UDI-DI: 611800222ADRUPVZ
Luer 3 parts Syringe 1 ml, 3 ml, 5 ml, 10 ml, 20 ml	18G 21G 22G	
	23G 24G 25G	
	28G 29G 30G	

Product Components

Name	Raw material
Barrel	PP - Polypropylene
Plunger	PP - Polypropylene
Gasket	Thermoplastic Elastomer, Latex free
Hypodermic Needle	Stainless Steel, SS304, PP – Polypropylene
Packaging	Medical grade paper and transparent film for EtO (individually blister packed)

Product Description

- A sterile device consisting of a calibrated barrel (cylinder) with plunger (attached with a gasket allowing smooth plunger movement) intended for the injection immediately after filling, a hypodermic needle is included, either attached or detached (capped for user protection).
- The syringe is typically made of plastic and Thermoplastic Elastomer materials and provided with re-use prevention feature which automatically locks the plunger upon use.
- This device is for single use only.

Measuring Function

This syringe includes a calibrated scale on the barrel to support accurate dose preparation. The graduations are expressed in millilitres (mL), in line with legal units under Directive 80/181/EEC.

The measuring function is validated according to ISO 7886-1. The graduation scale provides an accuracy within $\pm 5\%$ of the indicated volume for syringes ≥ 1 mL. Accuracy, precision and stability are maintained for the intended single-use duration when the syringe is intact, sterile and used as described in this IFU.

To achieve correct dosing, the user should read the scale at eye level, remove air bubbles before setting the final dose, and avoid any forward movement of the plunger before filling, as this activates the auto-disable mechanism.

Accuracy and Performance

The scale is designed and verified to provide consistent measurement accuracy for the intended purpose of preparing medication for immediate injection. The accuracy limits are established using validated test methods and are maintained throughout the intended single use.

HYPODERMIC AD/RUP SYRINGE (WITH/WITHOUT NEEDLE)

User Responsibilities

To achieve accurate dosing, the user must:

- Read the graduation marks at eye level.
- Remove air bubbles before setting the final dose.
- Avoid pushing the plunger forward before medication filling, as this will activate the auto-disable mechanism.
- Not use the syringe for blood withdrawal, fluid aspiration beyond immediate medication filling, or any purpose outside the stated intended use.

Limitations of the Measuring Function

The device is intended for one-time dose preparation immediately before injection.

It is not designed for long-term liquid storage or repeated measurement.

The accuracy of the scale is valid only under normal room-temperature use and when the syringe is intact, undamaged, and sterile.

Intended Use

Hypodermic AD/RUP Syringe is intended for the injection immediately after filling, it is mostly used for the administration of medications.

Intended User

The medical device shall be used by Healthcare Professional.

Patient Population

Hypodermic AD/RUP Syringe is designed for use with individuals of all ages, from newborns to adults and the elderly.


i Instructions for Use

DON'T PUSH THE PLUNGER FORWARD BEFORE FILLING THE MEDICATION.

Step 1: Open the pack and take out the syringe.

Do not push the plunger forward. Otherwise, the syringe will get locked.

Step 2: Insert the needle into the vial and withdraw the medication.


Step 3: Flick the barrel to move any air bubbles to the top.


Push the plunger to the desired dosage and expel the air bubbles.

HYPODERMIC AD/RUP SYRINGE (WITH/WITHOUT NEEDLE)


Step 4: Inject the medication into the patient until the full dose is delivered.

Step 5: Continue with the normal injection procedure as per clinical practice. Do not pull back the plunger after injection, as the AD mechanism will activate.

Step 6: For reconstitution, after Step 3, pull the plunger back fully, keeping the needle pointing downward to draw air into the syringe.

Step 7: Release the total fluid into the injection vial, keeping the vial downward. Then, withdraw the medication as in Steps 4, 3, 2, and 5.

Note: The plunger must not cross the ring, or the AD (Auto Disable) mechanism will activate.

Step 8: Dispose of the syringe in a sharps container.

The electronic version of "IFU/ADRUP/001 - Instructions for use" is available at www.immedicare.ma.

Compatibility with other medical devices

Suitable for use with connectors "EN ISO 80369-7 Small-bore connectors for liquids and gases in healthcare applications - Part 7: Connectors for intravenous or hypodermic applications".

Indications

- Hypodermic AD/RUP Syringes are devices used for administering medications.
- Syringes are designed to be used for intravenous, intramuscular, subcutaneous.

Contraindications

- Hypodermic AD/RUP Syringe should not be used in more than one patient. The product can be used in a one-time sterile condition with a single patient.
- Cleaning or re-sterilization is not possible. It is not intended for use other than its stated intended use.
- Dispose of as medical waste after use.

Cautions

- Do not use the product if the package is opened and damaged.

HYPODERMIC AD/RUP SYRINGE (WITH/WITHOUT NEEDLE)

- Visually inspect that the content of the package is complete, and the package is undamaged.
- Do not use the product after the expiration date printed on the primary packaging.
- After use, dispose of the product as medical waste.
- The product is for single use only. It cannot be reused or re-sterilized.
- Do not touch the needle tip to prevent contamination/needle sticks. Skin puncture with a contaminated needle may result in serious illnesses such as hepatitis, AIDS and known/unknown diseases.
- If the sterile packaging is not used immediately after opening, there is a risk of contamination. Product sterility cannot be maintained, do not use the product.
- Reuse of the product may cause infections, cross-contamination and sepsis.
- Sterile syringes are designed to be used immediately after being filled, so do not store the medicine for a long time.
- If the wrapping or cap on the end of the syringe is broken or missing, discard it in your sharps container and use a new needle and syringe.

Disclaimer

- The medical device is supplied with limited functionality. I'M MEDICARE does not warrant any foreseeable risks resulting from improper use or misuse of the device.
- The device should only be used according to its instructions for use. If the device is used outside the instructions described, the manufacturer or supplier must be informed of any off-label use.
- If any serious incident occurs, it should be reported to the manufacturer and the competent authority.

For any feedback or inquiries, customers/users can contact us at:

Email: contact@immedicare.ma

Customer Care No: +212 524 32 08 77

Adverse Event

- The device is supplied with limited functionality, and there are no foreseeable risks associated with its use.
- Improper use may cause redness, swelling, bruising, or other adverse effects.

Device Malfunction/Failure

In case of a device malfunction or failure, follow these steps:

- Discontinue use immediately if you notice any leakage, difficulty in activation, or failure to perform as intended.
- Do not attempt to repair the device yourself.
- Dispose of the device safely following the provided disposal guidelines.
- Report the malfunction to the manufacturer or an authorized representative for further investigation.

Storage Conditions

- Store in a cool and dry place.
- Keep it away from direct sunlight and moist environment.
- The shelf life of product is 5 years.

HYPODERMIC AD/RUP SYRINGE (WITH/WITHOUT NEEDLE)

Description of Symbols Used

Symbol	Description	Symbol	Description
	Manufactured by		Keep away from sunlight
	Date of manufacturing		Keep dry
	Use-by / Expiry date		Temperature limitation
LOT	Lot number / Batch number		Moisture limitation
REF	Product reference		Do not use if the package is damaged
MD	Medical device		Single sterile barrier system
UDI	Unique device identifier		Do not reuse
EC REP	European Authorized Representative		Latex free
CE	CE marking		Non-pyrogenic
STERILE EO	Sterilized by Ethylene oxide		Consult for IFU
	Do not re-sterilize		Do not use hook
	Caution / Warning (Read IFU before use)		Stock and store in upright condition
	Fragile. Handle with care		Re-use prevention

I'M MEDICARE SARL

Address: Parc Industriel Sidi Bou Othmane, Lot 21 Commune Sidi Bou Othmane, Benguerir, Marrakech, 43303, Morocco

Email: contact@immedicare.ma

Phone number: +212 631 700 900

Fax: +212 521 198 808

Website: www.immedicare.ma

EC REP

Name: EUCEREP B.V.

Address: Roald Dahllaan 33, 5629MC, Eindhoven, The Netherlands

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Nome: EUCEREP B.V.