



E008: Injection devices for immunization

Product name:	0.5ML AD SYRINGE FOR FIXED DOSE IMMUNIZATION
IMD-PQS Code:	E008-125
Status:	Prequalified
Date of acceptance:	20 Aug 2025
Product description:	0.5ML AD SYRINGE FOR FIXED DOSE IMMUNIZATION
Manufacturer:	IM Medicare
Manufacturers Reference:	0.5ML AD SYRINGE FOR FIXED DOSE IMMUNIZATION
Country of Manufacture:	Morocco
Address:	PARC INDUSTRIEL SIDI BOU OTHMANE, LOT 21 COMMUNE SIDI BOU OTHMANE, BENGUERIR Marrakech-Safi Morocco
Telephone:	+212 631 700 900
Email:	contact@immedicare.ma
Website address:	www.immedicare.ma
Valid until:	31 May 2026

Specifications

PRODUCT SPECIFICATIONS - MAIN

Vaccine capacity (mL)	0.50	Graduations (ml)	0.50
Markings	Under review	Syringe material(s)	Other
Syringe Material(s) (Other)	POLYPROPYLENE, THERMOPLASTIC ELASTOMER	Packaging & Sterility Assurance	Yes
Other needle size options	Yes	AD Number pieces	3
Auto-disabling at:	Start of injection	Needle Size	23Gx1" (0.60 x 25 mm)
AD Primary packaging	Blister pack		

Comments	MEDICATION FIXED DOSE IMMUNIZATION AD Syringe with needle is used for immunization with a fixed dose. It is used for administering medications /vaccines. It consists of a transparent barrel to ensure controlled administration of medicines with graduation, a plunger with seal/gasket and a fixed hypodermic needle (capped for user protection).
-----------------	--

PRICE AND SHIPPING

Shipping volume (m3)	0.057	Weight per Unit (Kg)	6.7000
Pieces per package	1800	Minimum order (unit)	180000.00
Price 1	On request	Year of base price	2025
Incoterms	FOB		

QUALITY STANDARD

Quality Standard	ISO 13485
-------------------------	-----------

Current PQS status:	<i>Prequalified</i>	Valid until:	<i>31 May 2026</i>
----------------------------	---------------------	---------------------	--------------------

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased.